

## Product datasheet for **TA315526**

### CSGLCAT (CHPF2) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human CSGLCAT.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	chondroitin polymerizing factor 2
Database Link:	<a href="#">NP_061888</a> <a href="#">Entrez Gene 54480 Human</a> <a href="#">Q9P2E5</a>
Synonyms:	chPF-2; ChSy-3; CSGLCA-T; CSGLcAT
Note:	CSGLCAT antibody detects endogenous levels of total CSGLCAT protein.
Protein Pathways:	Chondroitin sulfate biosynthesis, Metabolic pathways



[View online »](#)

## Product images:



Western blot analysis of extracts from Jurkat cells, using CSGLCAT antibody. The lane on the right is treated with the synthesized peptide.