

Product datasheet for **TA315785**

NOX5 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from C-terminal of human NOX5.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	NADPH oxidase, EF-hand calcium binding domain 5
Database Link:	NP_001171708 Entrez Gene 79400 Human Q96PH1
Synonyms:	MGC149776; MGC149777
Note:	NOX5 antibody detects endogenous levels of total NOX5 protein.
Protein Families:	Ion Channels: Other, Transmembrane



[View online »](#)

Product images:



Western blot analysis of extracts from Jurkat cells, using NOX5 antibody. The lane on the right is treated with the synthesized peptide.