

Product datasheet for **TA315903**

ANKRD20A1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human A20A1.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ankyrin repeat domain 20 family member A1
Database Link:	NP_115626 Entrez Gene 84210 Human Q5TYW2
Synonyms:	ANKRD20A
Note:	A20A1 antibody detects endogenous levels of total A20A1 protein.



[View online »](#)

Product images:



Western blot analysis of extracts from 3T3 cells, using A20A1 antibody. The lane on the right is treated with the synthesized peptide.