

## Product datasheet for **TA315939**


### ABHD14B Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 1:50~1:100, ELISA: 1:5000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human ABHD14B.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	abhydrolase domain containing 14B
Database Link:	<a href="#">NP_116139</a> <a href="#">Entrez Gene 76491 Mouse</a> <a href="#">Entrez Gene 300983 Rat</a> <a href="#">Entrez Gene 84836 Human</a> <a href="#">Q96IU4</a>
Synonyms:	CIB; HEL-S-299
Note:	ABHD14B antibody detects endogenous levels of total ABHD14B protein.



[View online »](#)

**Product images:**

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ABHD14B antibody. The picture on the right is treated with the synthesized peptide.