

Product datasheet for **TA315968**

ATPBD3 (CTU1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:20000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from C-terminal of human ATPB3.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	cytosolic thiouridylase subunit 1
Database Link:	NP_660275 Entrez Gene 90353 Human Q7Z7A3
Synonyms:	ATPBD3; NCS6
Note:	ATBP3 antibody detects endogenous levels of total ATPB3 protein.



[View online »](#)

Product images:



Western blot analysis of extracts from LOVO cells, using ATPBP3 antibody. The lane on the right is treated with the synthesized peptide.