


Product datasheet for **TA315969**

GADD45GIP1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 1:50~1:100, ELISA: 1:40000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human GADD45GIP1.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	GADD45G interacting protein 1
Database Link:	NP_443082 Entrez Gene 90480 Human Q8TAE8
Synonyms:	CKBBP2; CKbetaBP2; CRIF1; MRP-L59; PLINP; PLINP-1; Plinp1; PRG6
Note:	GADD45GIP1 antibody detects endogenous levels of total GADD45GIP1 protein.

[View online »](#)

Product images:

Immunohistochemistry analysis of paraffin-embedded human brain tissue using GADD45GIP1 antibody. The picture on the right is treated with the synthesized peptide.