

Product datasheet for **TA315999**


MUC16 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF
Recommended Dilution:	IF: 1:100~1:500, ELISA: 1:40000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human MUC16.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mucin 16, cell surface associated
Database Link:	NP_078966 Entrez Gene 94025 Human Q8WXI7
Synonyms:	CA125
Note:	MUC16 antibody detects endogenous levels of total MUC16 protein.



[View online »](#)

Product images:

Immunofluorescence analysis of HepG2 cells, using MUC16 antibody. The picture on the right is treated with the synthesized peptide.