

## Product datasheet for **TA316067**

### SLC24A4 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human SLC24A4.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	solute carrier family 24 member 4
Database Link:	<a href="#">NP_705932</a> <a href="#">Entrez Gene 238384 MouseEntrez Gene 314396 RatEntrez Gene 123041 Human Q8NFF2</a>
Synonyms:	AI2A5; NCKX4; SHEP6; SLC24A2
Note:	SLC24A4 antibody detects endogenous levels of total SLC24A4 protein.
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

## Product images:



Western blot analysis of extracts from RAW264.7 cells, using SLC24A4 antibody. The lane on the right is treated with the synthesized peptide.