

## Product datasheet for **TA316469**


### OR5B3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF
Recommended Dilution:	IF: 1:100~1:500, ELISA: 1:10000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human OR5B3.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	olfactory receptor family 5 subfamily B member 3
Database Link:	<a href="#">NP_001005469</a> <a href="#">Entrez Gene 441608 Human</a> <a href="#">Q8NH48</a>
Synonyms:	OR5B13; OR11-239; OST129
Note:	OR5B3 antibody detects endogenous levels of total OR5B3 protein.
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction



[View online »](#)

**Product images:**

Immunofluorescence analysis of COS-7 cells, using OR5B3 antibody. The picture on the right is treated with the synthesized peptide.