

Product datasheet for **TA316673**


DPEP1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (15 µg/ml)
Reactivity:	Human (Predicted: Monkey)
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Renal Dipeptidase / DPEP1 antibody was raised against synthetic 14 amino acid peptide from internal region of human DPEP1. Percent identity with other species by BLAST analysis: Human (100%); Gorilla, Gibbon, Monkey (93%); Panda, Rabbit, Horse (86%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	dipeptidase 1 (renal)
Database Link:	NP_004404 Entrez Gene 699725 Monkey Entrez Gene 1800 Human P16444
Synonyms:	MBD1; MDP; RDP
Note:	Specific for Human DPEP1. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Protein Families:	Protease



[View online »](#)

Product images:

Anti-DPEP1 antibody IHC of human colon. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.



Anti-DPEP1 antibody IHC of human uterus. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.