

Product datasheet for **TA316700**


Embigin homolog (EMB) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (10 µg/ml)
Reactivity:	Chimpanzee, Gorilla, Human, Monkey
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	EMB/ Embigin antibody was raised against synthetic 17 amino acid peptide from C-Terminus of human EMB. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Monkey (100%); Orangutan (94%); Gibbon, Marmoset, Bat (88%); Bovine, Dog (82%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	embigin
Database Link:	NP_940851 Entrez Gene 702068 Monkey Entrez Gene 133418 Human Q6PCB8
Synonyms:	GP70
Note:	Specific for Human EMB. BLAST analysis of the peptide immunogen showed no homology with other human proteins, except RALGPS1 (53%).
Protein Families:	Transmembrane



[View online »](#)

Product images:

Anti-EMB antibody IHC of human placenta. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.



Anti-EMB antibody IHC of human stomach. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.