

## **Product datasheet for TA317155**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## PDE8A Rabbit Polyclonal Antibody

**Product data:** 

**Product Type:** Primary Antibodies

Applications: IHC

**Recommended Dilution:** IHC-P (10 µg/ml)

**Reactivity:** Dog, Gorilla, Human, Pig, Rabbit, Gibbon, Horse (Predicted: Monkey)

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**Immunogen:** PDE8A antibody was raised against synthetic 16 amino acid peptide from near N-terminus of

human PDE8A. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon,

Monkey, Dog, Panda, Horse, Rabbit, Pig (100%); Marmoset (94%); Mouse, Rat, Bovine,

Elephant, Opossum (88%); Platypus, Zebrafish (81%).

**Formulation:** PBS, 0.1% sodium azide.

**Concentration:** lot specific

**Purification:** Immunoaffinity Chromatography

**Conjugation:** Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** phosphodiesterase 8A

Database Link: NP 775659

Entrez Gene 488752 DogEntrez Gene 694750 MonkeyEntrez Gene 5151 Human

O60658

Synonyms: 5'-cyclic phosphodiesterase 8A; cAMP-specific cyclic nucleotide phosphodiesterase 8A;

FLJ16150; high-affinity cAMP-specific and IBMX-insensitive 3'; HsT19550;

OTTHUMP00000192898; phosphodiesterase 8A

Note: Specific for Human PDE8A. BLAST analysis of the peptide immunogen showed no homology

with other human proteins.

**Protein Families:** Druggable Genome






**Protein Pathways:** 

Progesterone-mediated oocyte maturation, Purine metabolism

## **Product images:**



Anti-PDE8A antibody IHC of human prostate. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.