

Product datasheet for **TA317241**


QSOX2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (20 µg/ml)
Reactivity:	Gorilla, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	QSOX2 antibody was raised against synthetic 15 amino acid peptide from internal region of human QSOX2. Percent identity with other species by BLAST analysis: Human, Gorilla (100%); Marmoset (87%); Mouse, Rat, Dog, Elephant, Horse (80%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	quiescin sulfhydryl oxidase 2
Database Link:	NP_859052 Entrez Gene 169714 Human Q6ZRP7
Synonyms:	QSCN6L1; SOXN
Note:	Specific for Human QSOX2. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Protein Families:	Transmembrane



[View online »](#)

Product images:

Anti-QSOX2 antibody IHC of human adrenal. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.



Anti-QSOX2 antibody IHC of human heart. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.