

Product datasheet for **TA317418**


Thrombopoietin (THPO) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	ELISA (1:000 - 1:1000), IHC-P (10 µg/ml)
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	THPO / TPO / Thrombopoietin antibody was raised against synthetic peptide from human THPO / Thrombopoietin.
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Protein G Column
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	thrombopoietin
Database Link:	NP_000451 Entrez Gene 7066 Human P40225
Synonyms:	MGDF; MKCSF; ML; MPLLG; THCYT1; TPO
Note:	Specific for Amino acids 35 to 50 of human THPO
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Hematopoietic cell lineage



[View online »](#)

Product images:

Anti-Thrombopoietin antibody IHC of human spleen. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 10 ug/ml.