

Product datasheet for **TA317770**

NALP2 (NLRP2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	IHC-P (5 µg/ml), WB (1 - 2 µg/ml)
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide from human NLRP2 / NALP2.
Formulation:	PBS, 0.02% sodium azide.
Concentration:	lot specific
Purification:	IgG Fractionation
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	NLR family, pyrin domain containing 2
Database Link:	NP_060322 Entrez Gene 232827 MouseEntrez Gene 684545 RatEntrez Gene 55655 Human Q9NX02
Synonyms:	CLR19.9; NALP2; NBS1; PAN1; PYPAF2
Note:	peptide corresponding to 12 amino acids near the carboxy-terminus of human NALP2



[View online »](#)