

## **Product datasheet for TA319422**

## **Asap1 Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type: Primary Antibodies** 

IF, IHC **Applications:** 

Recommended Dilution: ELISA: 1:4,000 - 1:16,000, WB: 1:500 - 1:2,000, IHC: 20-40 ug/ml

Reactivity: Mouse, Human, Chicken, Bovine, Dog, Rat, Chimpanzee

**Modifications:** Phospho-specific

Host: Rabbit Isotype: **IgG** 

Clonality: Polyclonal

Immunogen: This affinity purified antibody was prepared from whole rabbit serum produced by repeated

immunizations with a synthetic peptide corresponding to amino acids 777-787 of mouse

ASAP1 protein.

Formulation: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Concentration: lot specific

Unconjugated Conjugation:

Store at -20°C as received. Storage:

Stable for 12 months from date of receipt. Stability:

Gene Name: ArfGAP with SH3 domain, ankyrin repeat and PH domain1

Database Link: NP 034156

Entrez Gene 50807 HumanEntrez Gene 314961 RatEntrez Gene 482044 DogEntrez Gene

13196 Mouse

Q9QWY8

Synonyms: AMAP1; CENTB4; DDEF1; DEF-1; KIAA1249; PAG2; PAP; ZG14P



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Note:

This antibody is suitable for Cancer, Immunology and Nuclear Signaling research. ASAP1 (also known as AMAP1, 130-kDa phosphatidylinositol 4,5-biphosphate-dependent ARF1 GTPase-activating protein, PIP2-dependent ARF1 GAP,ADP-ribosylation factor-directed GTPase-activating protein 1, ARF GTPase-activating protein 1, Development and differentiation-enhancing factor 1, Differentiation-enhancing factor 1, DEF-1) is an Arf-directed GTPase activating protein that is a substrate for the kinases Src and FAK and has been implicated in the regulation of membrane traffic, focal adhesions and invadopodia/podosomes. Phosphorylation of ASAP1 at tyrosine 782 has been found to affect enzymatic and some biological activities, including the function of invadopodia. ASAP1 is expressed in many tissues but is most abundant in the testis, brain, lung and spleen. A heightened expression was seen in the adipose tissue from obese (ob) and diabetic (db) animals. Multiple transcript variants have been reported for this protein.

## **Product images:**





Anti-ASAP1 pY782 antibody was used at 20 ?g/ml to detect signal in a variety of tissues including multi-human, multi-brain and multi-cancer slides. This image shows moderate intracellular positive staining in epidermal keratinocytes in human skin at 40X. Tissue was formalin-fixed and paraffin embedded. The image shows localization of the antibody as the precipitated red signal, with a hematoxylin purple nuclear counterstain. Personal Communi-cation, Tina Roush, LifeSpanBiosciences, Seattle, WA.

Immunofluorescent microscopy using Affinity Purified anti-ASAP1 pY782 antibody shows detection of phosphorylated ASAP1 present in mouse NIH3T3 cells transfected with activated Src. Specific staining is not present when antibody is pre-incubated with the immunizing peptide prior to reaction with cells. Personal Communication. Paul Randazzo, NIH, CCR, Bethesda, MD.