

Product datasheet for **TA321073**

Glucokinase (GCK) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 3-17 amino acids of Human glucokinase (hexokinase 4)
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	glucokinase
Database Link:	NP_000153 Entrez Gene 24385 Rat Entrez Gene 103988 Mouse Entrez Gene 2645 Human P35557




[View online »](#)

- Background:** Hexokinases phosphorylate glucose to produce glucose-6-phosphate; the first step in most glucose metabolism pathways. Alternative splicing of this gene results in three tissue-specific forms of glucokinase; one found in pancreatic islet beta cells and two found in liver. The protein localizes to the outer membrane of mitochondria. In contrast to other forms of hexokinase; this enzyme is not inhibited by its product glucose-6-phosphate but remains active while glucose is abundant. Mutations in this gene have been associated with non-insulin dependent diabetes mellitus (NIDDM); maturity onset diabetes of the young; type 2 (MODY2) and persistent hyperinsulinemic hypoglycemia of infancy (PHHI).
- Synonyms:** FGQTL3; GK; GLK; HHF3; HK4; HKIV; HXKP; LGLK; MODY2
- Protein Families:** Druggable Genome
- Protein Pathways:** Amino sugar and nucleotide sugar metabolism, Galactose metabolism, Glycolysis / Gluconeogenesis, Insulin signaling pathway, Maturity onset diabetes of the young, Metabolic pathways, Starch and sucrose metabolism, Type II diabetes mellitus

Product images:



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA321073 (GCK Antibody) at dilution 1/40 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA321073 (GCK Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)