

Product datasheet for TA321080

PDX1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 200-1000

WB positive control: Hela and Jurkat cells

IHC: 15-50

Positive control: Human gasrtic cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Synthetic peptide peptide corresponding to a region derived from 18-31 amino acids of

human pancreatic and duodenal homeobox 1

Formulation: PBS pH7.3, 0.05% NaN3, 50% glycerol

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 31 kDa

Gene Name: pancreatic and duodenal homeobox 1

Database Link: NP 000200

Entrez Gene 18609 MouseEntrez Gene 29535 RatEntrez Gene 3651 Human

P52945



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background: The protein encoded by this gene is a transcriptional activator of several genes; including

insulin; somatostatin; glucokinase; islet amyloid polypeptide; and glucose transporter type 2. The encoded nuclear protein is involved in the early development of the pancreas and plays a major role in glucose-dependent regulation of insulin gene expression. Defects in this gene are a cause of pancreatic agenesis; which can lead to early-onset insulin-dependent diabetes

mellitus (NIDDM); as well as maturity onset diabetes of the young type 4 (MODY4).


Synonyms: GSF; IDX-1; IPF1; IUF1; MODY4; PAGEN1; PDX-1; STF-1

Protein Families: Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transcription

Factors

Protein Pathways: Maturity onset diabetes of the young, Type II diabetes mellitus

Product images:



Gel: 10%SDS-PAGE Lysate: 40 µg Lane 1-2: Hela cells Jurkat cells


Primary antibody: TA321080 (PDX1 Antibody) at

dilution 1/400

Secondary antibody: Goat anti rabbit IgG at


1/8000 dilution

Exposure time: 1 minute



Immunohistochemistry of paraffin-embedded Human gasrtic cancer tissue using TA321080 (PDX1 Antibody) at dilution 1/20 (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human gasrtic cancer tissue using TA321080 (PDX1 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: ×200)