

Product datasheet for **TA321144S**

Adenylate Kinase 1 (AK1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 200-1000 WB positive control: Human fetal muscle tissue
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to a region derived from 1-194 amino acids of Human Adenylate kinase isoenzyme 1
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	22 kDa
Gene Name:	adenylate kinase 1
Database Link:	NP_000467 Entrez Gene 11636 Mouse Entrez Gene 24183 Rat Entrez Gene 203 Human P00568
Background:	Adenylate kinase is an enzyme involved in regulating the adenine nucleotide composition within a cell by catalyzing the reversible transfer of phosphate group among adinine nucleotides. Three isozymes of adenylate kinase have been identified in vertebrates; adenylate isozyme 1 (AK1); 2 (AK2) and 3 (AK3). AK1 is found in the cytosol of skeletal muscle; brain and erythrocytes; whereas AK2 and AK3 are found in the mitochondria of other tissues including liver and heart. AK1 was identified because of its association with a rare genetic disorder causing nonspherocytic hemolytic anemia where a mutation in the AK1 gene was found to reduce the catalytic activity of the enzyme.




[View online »](#)

Synonyms: HTL-S-58j

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Purine metabolism

Product images:

Gel: 10%SDS-PAGE
Lysate: 40 µg
Lane: Human fetal muscle tissue
Primary antibody: [TA321144] (AK1 Antibody) at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 10 seconds