

Product datasheet for **TA321222**

HTR2B Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: 293T cells IHC: 25-100 Positive control: Human lung cancer Predicted cell location: Cytoplasm, Cell membrane
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to a region derived from 1-58 amino acids of Human 5-hydroxytryptamine receptor 2B
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54 kDa
Gene Name:	5-hydroxytryptamine receptor 2B
Database Link:	NP_000858 Entrez Gene 3357 Human P41595



[View online »](#)

Background:

This gene encodes one of the several different receptors for 5-hydroxytryptamine (serotonin) that belongs to the G-protein coupled receptor 1 family. Serotonin is a biogenic hormone that functions as a neurotransmitter; a hormone; and a mitogen. Serotonin receptors mediate many of the central and peripheral physiologic functions of serotonin; including regulation of cardiovascular functions and impulsive behavior. Population and family-based analyses of a minor allele (glutamine-to-stop substitution; designated Q20*) which blocks expression of this protein; and knockout studies in mice; suggest a role for this gene in impulsivity. However; other factors; such as elevated testosterone levels; may also be involved.

Synonyms:

5-HT(2B); 5-HT2B

Protein Families:

Druggable Genome, Transmembrane

Protein Pathways:

Calcium signaling pathway, Gap junction, Neuroactive ligand-receptor interaction

Product images:

Gel: 10%SDS-PAGE


Lysate: 20 µg

Lane: 293T cells


Primary antibody: TA321222 (HTR2B Antibody) at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA321222 (HTR2B Antibody) at dilution 1/20 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA321222 (HTR2B Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: ×200)