

Product datasheet for **TA321317**

CDC20 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human cervical cancer tissue IHC: 5-25 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to a region derived from 182-480 amino acids of Human Cell division cycle protein 20 homolog
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55 kDa
Gene Name:	cell division cycle 20
Database Link:	NP_001246 Entrez Gene 64515 Rat Entrez Gene 107995 Mouse Entrez Gene 991 Human Q12834
Background:	CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes; nuclear movement prior to anaphase and chromosome separation.
Synonyms:	bA276H19.3; CDC20A; p55CDC



[View online »](#)


Protein Families: Druggable Genome

Protein Pathways: Cell cycle, Oocyte meiosis, Ubiquitin mediated proteolysis


Product images:



Gel: 10%SDS-PAGE
Lysate: 30 μ g
Lane: Human cervical cancer tissue
Primary antibody: TA321317 (CDC20 Antibody) at dilution 1/600
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 1 minute



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA321317 (CDC20 Antibody) at dilution 1/10 (Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA321317 (CDC20 Antibody) at dilution 1/10, treated with fusion protein. (Original magnification: ×200)