

Product datasheet for **TA321324S**

Claudin 4 (CLDN4) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: HT-29 cells IHC: 25-100 Positive control: Human breast cancer Predicted cell location: Cytoplasm, Cell membrane
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 193-209 amino acids of Human claudin 4
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	22 kDa
Gene Name:	claudin 4
Database Link:	NP_001296 Entrez Gene 12740 Mouse Entrez Gene 1364 Human O14493



[View online »](#)

Background:

Claudin 4; also known as CLDN4; is a protein which in humans is encoded by the CLDN4 gene. It belongs to the group of claudins. This gene encodes an integral membrane protein; which belongs to the claudin family. The protein is a component of tight junction strands and may play a role in internal organ development and function during pre- and postnatal life. This gene is deleted in Williams-Beuren syndrome; a neurodevelopmental disorder affecting multiple systems.

Synonyms:


CPE-R; CPER; CPETR; CPETR1; hCPE-R; WBSCR8

Protein Families:


Druggable Genome, Transmembrane

Protein Pathways:


Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction

Product images:


Gel: 10+12%SDS-PAGE
 Lysate: 30 µg
 Lane: HT29 cells
 Primary antibody: [TA321324] (CLDN4 Antibody) at dilution 1/700
 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
 Exposure time: 30 seconds



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA321324] (CLDN4 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using [TA321324] (CLDN4 Antibody) at dilution 1/35, treated with synthetic peptide. (Original magnification: ×200)