

Product datasheet for **TA321391**

ARF6 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: HepG2 and 231 cells IHC: 10-50 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 162-175 amino acids of Human ADP-ribosylation factor 6
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	20 kDa
Gene Name:	ADP ribosylation factor 6
Database Link:	NP_001654 Entrez Gene 11845 Mouse Entrez Gene 79121 Rat Entrez Gene 382 Human P62330



[View online »](#)

Background:


This gene encodes a member of the human ARF gene family; which is part of the RAS superfamily. The ARF genes encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking and as activators of phospholipase D. The product of this gene is localized to the plasma membrane; and regulates vesicular trafficking; remodelling of membrane lipids; and signaling pathways that lead to actin remodeling. A pseudogene of this gene is located on chromosome 7.

Synonyms:


DKFZp564M0264

Protein Pathways:


Endocytosis, Fc gamma R-mediated phagocytosis

Product images:

Gel: 12%SDS-PAGE
Lysate: 40 µg
Lane 1-2: HepG2 cells
231 cells
Primary antibody: TA321391 (ARF6 Antibody) at dilution 1/300
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 1 minute



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA321391 (ARF6 Antibody) at dilution 1/15 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA321391 (ARF6 Antibody) at dilution 1/15, treated with synthetic peptide. (Original magnification: ×200)