

Product datasheet for **TA321493S**

NGF Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Hela cell lysate IHC: 50-100 Positive control: Human tonsil Predicted cell location: Secreted
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 80-94 amino acids of human nerve growth factor (beta polypeptide)
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	27 kDa
Gene Name:	nerve growth factor
Database Link:	NP_002497 Entrez Gene 18049 MouseEntrez Gene 310738 RatEntrez Gene 4803 Human P01138



[View online »](#)


Background: This gene is a member of the NGF-beta family and encodes a secreted protein which homodimerizes and is incorporated into a larger complex. This protein has nerve growth stimulating activity and the complex is involved in the regulation of growth and the differentiation of sympathetic and certain sensory neurons. Mutations in this gene have been associated with hereditary sensory and autonomic neuropathy; type 5 (HSAN5); and dysregulation of this gene's expression is associated with allergic rhinitis.

Synonyms: Beta-NGF; HSAN5; NGFB


Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Apoptosis, MAPK signaling pathway, Neurotrophin signaling pathway


Product images:



Gel: 12%SDS-PAGE
Lysate: 40 µg
Lane: HeLa cell lysate
Primary antibody: [TA321493] (NGF Antibody) at dilution 1/700
Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution
Exposure time: 3 minutes



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA321493] (NGF Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA321493] (NGF Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: $\times 200$)