

## **Product datasheet for TA321538**

## **SHP2 (PTPN11) Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 500-2000

WB positive control: NIH/3T3, 293T, 231 cell lysates

IHC: 50-200

Positive control: Human esophagus cancer

Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**Immunogen:** Fusion protein corresponding to C terminal 250 amino acids of human protein tyrosine

phosphatase, non-receptor type 11

**Formulation:** PBS pH7.3, 0.05% NaN3, 50% glycerol

**Concentration:** lot specific

**Purification:** Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 68 kDa

**Gene Name:** protein tyrosine phosphatase, non-receptor type 11

Database Link: NP 002825

Entrez Gene 19247 MouseEntrez Gene 25622 RatEntrez Gene 5781 Human

Q06124



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



## Background:

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth; differentiation; mitotic cycle; and oncogenic transformation. This PTP contains two tandem Src homology-2 domains; which function as phospho-tyrosine binding domains and mediate the interaction of this PTP with its substrates. This PTP is widely expressed in most tissues and plays a regulatory role in various cell signaling events that are important for a diversity of cell functions; such as mitogenic activation; metabolic control; transcription regulation; and cell migration. Mutations in this gene are a cause of Noonan syndrome as well as acute myeloid leukemia. Two transcript variants encoding different isoforms have been found for this gene.

Synonyms: BPTP3; CFC; JMML; METCDS; NS1; PTP-1D; PTP2C; SH-PTP2; SH-PTP3; SHP2

**Protein Families:** Druggable Genome, Phosphatase


**Protein Pathways:** Adipocytokine signaling pathway, Chronic myeloid leukemia, Epithelial cell signaling in

Helicobacter pylori infection, Jak-STAT signaling pathway, Leukocyte transendothelial

migration, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Renal cell

carcinoma

## **Product images:**



Gel: 8%SDS-PAGE Lysate: 40 µg Lane 1-3: NIH/3T3 293T

231 cell lysates

Primary antibody: TA321538 (PTPN11 Antibody)


at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at


1/5000 dilution

Exposure time: 20 seconds






Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA321538 (PTPN11 Antibody) at dilution 1/60 (Original magnification: ×200)




Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA321538 (PTPN11 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using TA321538 (PTPN11 Antibody) at dilution 1/60 (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using TA321538 (PTPN11 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: ×200)