

Product datasheet for TA321767

PTP epsilon (PTPRE) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 15-50

Positive control: Human ovarian cancer

Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Fusion protein corresponding to a region derived from 426-689 amino acids of human

protein tyrosine phosphatase, receptor type, E

Formulation: PBS pH7.3, 0.05% NaN3, 50% glycerol

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: protein tyrosine phosphatase, receptor type E

Database Link: NP 006495

Entrez Gene 19267 MouseEntrez Gene 114767 RatEntrez Gene 5791 Human

P23469



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Background:

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth; differentiation; mitotic cycle; and oncogenic transformation. Two alternatively spliced transcript variants of this gene have been reported; one of which encodes a receptor-type PTP that possesses a short extracellular domain; a single transmembrane region; and two tandem intracytoplasmic catalytic domains; Another one encodes a PTP that contains a distinct hydrophilic N-terminus; and thus represents a nonreceptor-type isoform of this PTP. Studies of the similar gene in mice suggested the regulatory roles of this PTP in RAS related signal transduction pathways; cytokines induced SATA signaling; as well as the activation of voltage-gated K+ channels.


Synonyms: HPTPE; PTPE; R-PTP-EPSILON

Protein Families: Druggable Genome, Phosphatase, Transmembrane

Product images:




Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA321767 (PTPRE Antibody) at dilution 1/10 (Original magnification: ×200)




Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA321767 (PTPRE Antibody) at dilution 1/10, treated with fusion protein. (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA321767 (PTPRE Antibody) at dilution 1/10 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA321767 (PTPRE Antibody) at dilution 1/10, treated with fusion protein. (Original magnification: ×200)