

Product datasheet for **TA321846**

Adiponectin Receptor 1 (ADIPOR1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 200-1000 WB positive control: Mouse heart tissue
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to a region derived from 1-136 amino acids of human adiponectin receptor 1
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	43 kDa
Gene Name:	adiponectin receptor 1
Database Link:	NP_057083 Entrez Gene 72674 Mouse Entrez Gene 51094 Human Q96A54
Background:	This gene encodes a protein which acts as a receptor for adiponectin; a hormone secreted by adipocytes which regulates fatty acid catabolism and glucose levels. Binding of adiponectin to the encoded protein results in activation of an AMP-activated kinase signaling pathway which affects levels of fatty acid oxidation and insulin sensitivity. A pseudogene of this gene is located on chromosome 14.?
Synonyms:	ACDCR1; CGI-45; CGI45; PAQR1; TESBP1A




[View online »](#)

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Adipocytokine signaling pathway

Product images:



Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane: Mouse heart tissue
Primary antibody: TA321846 (ADIPOR1 Antibody)
at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 30 seconds