

Product datasheet for **TA321876S**


RAB25 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-150 Positive control: Human liver cancer Predicted cell location: Cell membrane, Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 127-140 amino acids of Human Ras-related protein Rab-25
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	RAB25, member RAS oncogene family
Database Link:	NP_065120 Entrez Gene 53868 Mouse Entrez Gene 57111 Human P57735
Background:	RAB proteins; such as RAB25; are members of the RAS superfamily of small GTPases that are involved in membrane trafficking. Members of the RAB11 subfamily; including RAB25; control the return of internalized membrane-associated moieties to the cell surface.
Synonyms:	CATX-8; RAB11C
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA321876] (RAB25 Antibody) at dilution 1/50. (Original magnification: $\times 200$)