

Product datasheet for **TA322091S**

ACY3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Mouse kidney tissue IHC: 25-100 Positive control: Human cervical cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	35 kDa
Gene Name:	aminoacylase 3
Database Link:	NP_542389 Entrez Gene 71670 Mouse Entrez Gene 91703 Human Q96HD9

Background: Acy3, which encodes a deduced 318-amino acid protein with a calculated molecular mass of about 35.2 kD. Acy3 has no membrane domain, but it contains consensus sites for N-glycosylation, tyrosine sulfation, tyrosine and threonine phosphorylation, and myristoylation. Northern blot analysis detected high expression of an approximately 1.4-kb Acy3 transcript in kidney, followed by liver. Weaker expression was detected in heart, small intestine, brain, lung, testis, and stomach. Testis expressed an additional transcript of about 2.2 kb.




[View online »](#)


Synonyms: ACY-3; HCBP1

Protein Pathways: Alanine, aspartate and glutamate metabolism, Histidine metabolism


Product images:




Gel: 10%SDS-PAGE
Lysate: 40 µg
Lane: Mouse kidney tissue
Primary antibody: [TA322091] (ACY3 Antibody) at dilution 1/500
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 5 seconds




Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA322091] (ACY3 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA322091] (ACY3 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA322091] (ACY3 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA322091] (ACY3 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: ×200)