

## Product datasheet for **TA322132S**

### alpha Internexin (INA) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1000-5000 WB positive control: Mouse brain tissue IHC: 50-200 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 486-499 amino acids of Human internexin neuronal intermediate filament protein, alpha
Formulation:	PBS pH7.3, 0.05% NaN <sub>3</sub> , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55 kDa
Gene Name:	internexin neuronal intermediate filament protein alpha
Database Link:	<a href="#">NP_116116</a> <a href="#">Entrez Gene 226180 Mouse</a> <a href="#">Entrez Gene 9118 Human</a> <a href="#">Q16352</a>
Background:	Neurofilaments are type IV intermediate filament heteropolymers composed of light; medium; and heavy chains. Neurofilaments comprise the axoskeleton and they functionally maintain the neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. This gene is a member of the intermediate filament family and is involved in the morphogenesis of neurons.?




[View online »](#)

Synonyms: NEF5; NF-66; TXBP-1


### Product images:



Gel: 10%SDS-PAGE  
Lysate: 50 µg  
Lane: Mouse brain tissue  
Primary antibody: [TA322132] (INA Antibody) at dilution 1/800  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
Exposure time: 2 seconds



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA322132] (INA Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA322132] (INA Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: x200)