

## Product datasheet for **TA322233**

### Adenylate cyclase 1 (ADCY1) Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 100-300 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 1057-1070 amino acids of human adenylate cyclase 1 (brain)
Formulation:	PBS pH7.3, 0.05% NaN <sub>3</sub> , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	adenylate cyclase 1 (brain)
Database Link:	<a href="#">NP_066939</a> <a href="#">Entrez Gene 107 Human</a> <a href="#">Q08828</a>
Background:	This gene encodes a form of adenylate cyclase expressed in brain. A similar protein in mouse is involved in pattern formation of the brain.
Synonyms:	AC1; DFNB44
Protein Families:	Druggable Genome, Transmembrane




[View online »](#)

**Protein Pathways:**


Calcium signaling pathway, Chemokine signaling pathway, Dilated cardiomyopathy, Gap junction, GnRH signaling pathway, Long-term potentiation, Melanogenesis, Oocyte meiosis, Progesterone-mediated oocyte maturation, Purine metabolism, Vascular smooth muscle contraction

**Product images:**


Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA322233 (ADCY1 Antibody) at dilution 1/60 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA322233 (ADCY1 Antibody) at dilution 1/60, treated with synthetic peptide. (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA322233 (ADCY1 Antibody) at dilution 1/60 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA322233 (ADCY1 Antibody) at dilution 1/60, treated with synthetic peptide. (Original magnification:  $\times 200$ )