

## **Product datasheet for TA322241**

## **SCN2A Rabbit Polyclonal Antibody**

## **Product data:**

**Product Type:** Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 25-100

Positive control: Human thyroid cancer

Predicted cell location: Cytoplasm

Reactivity: Human, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**Immunogen:** Synthetic peptide corresponding to a region derived from 30-43 amino acids of human

sodium channel, voltage-gated, type II, alpha subunit

Formulation: PBS pH7.3, 0.05% NaN3, 50% glycerol

**Concentration:** lot specific

**Purification:** Antigen affinity purification

**Conjugation:** Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** sodium voltage-gated channel alpha subunit 2

Database Link: NP 066287

Entrez Gene 24766 RatEntrez Gene 6326 Human

Q99250



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background:

Voltage-gated sodium channels are transmembrane glycoprotein complexes composed of a large alpha subunit with 24 transmembrane domains and one or more regulatory beta subunits. They are responsible for the generation and propagation of action potentials in neurons and muscle. This gene encodes one member of the sodium channel alpha subunit gene family. It is heterogeneously expressed in the brain; and mutations in this gene have been linked to several seizure disorders. Several alternatively spliced transcript variants of this gene have been described; but the full-length nature of some of these variants has not been determined.


Synonyms:

BFIC3; BFIS3; BFNIS; EIEE11; HBA; HBSCI; HBSCII; Na(v)1.2; NAC2; Nav1.2; SCN2A1; SCN2A2


**Protein Families:** 

Druggable Genome, Ion Channels: Sodium, Transmembrane

## **Product images:**




Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA322241 (SCN2A Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA322241 (SCN2A Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA322241 (SCN2A Antibody) at dilution 1/25 (Original magnification: ×200)

Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using TA322241 (SCN2A Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: ×200)