

Product datasheet for **TA322301S**

ACSM5 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: 293T cell IHC: 25-100 Positive control: Human thyroid cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	acyl-CoA synthetase medium-chain family member 5
Database Link:	NP_060358 Entrez Gene 272428 Mouse Entrez Gene 54988 Human Q6NUN0
Background:	Fatty acids are incorporated into membranes and signaling molecules and have roles in energy storage and metabolism. These essential functions require activation of the fatty acid by acyl-coenzyme A (CoA) synthetases, such as ACSM5, which form an activating thioester linkage between the fatty acid and CoA. Has medium-chain fatty acid:CoA ligase activity with broad substrate specificity (in vitro). Acts on acids from C4 to C11 and on the corresponding 3-hydroxy- and 2,3- or 3,4-unsaturated acids (in vitro).
Synonyms:	FLJ20581




[View online »](#)

Protein Pathways: Butanoate metabolism, Metabolic pathways


Product images:




Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane: 293T cell
Primary antibody: [TA322301] (ACSM5 Antibody) at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 30 seconds




Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA322301] (ACSM5 Antibody) at dilution 1/35 (Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA322301] (ACSM5 Antibody) at dilution 1/35, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using [TA322301] (ACSM5 Antibody) at dilution 1/35 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using [TA322301] (ACSM5 Antibody) at dilution 1/35, treated with fusion protein. (Original magnification: $\times 200$)