

## Product datasheet for **TA322338**


### A4GNT Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to a region derived from 28-340 amino acids of human alpha-1,4-N-acetylglucosaminyltransferase
Formulation:	PBS pH7.3, 0.05% NaN <sub>3</sub> , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	alpha-1,4-N-acetylglucosaminyltransferase
Database Link:	<a href="#">NP_057245</a> <a href="#">Entrez Gene 51146 Human</a> <a href="#">Q9UNA3</a>
Background:	This gene encodes a protein from the glycosyltransferase 32 family. The enzyme catalyzes the transfer of N-acetylglucosamine (GlcNAc) to core 2 branched O-glycans. It forms a unique glycan; GlcNAc $\alpha$ 1 $\rightarrow$ 4Gal $\beta$ $\rightarrow$ R and is largely associated with the Golgi apparatus membrane.
Synonyms:	alpha4GnT
Protein Families:	Transmembrane




[View online »](#)

**Product images:**


Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA322338 (A4GNT Antibody) at dilution 1/25 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA322338 (A4GNT Antibody) at dilution 1/25, treated with fusion protein. (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA322338 (A4GNT Antibody) at dilution 1/25 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using TA322338 (A4GNT Antibody) at dilution 1/25, treated with fusion protein. (Original magnification:  $\times 200$ )