

Product datasheet for **TA322417S**

Claudin 14 (CLDN14) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 200-1000 WB positive control: Mouse pancreas tissue
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 227-239 amino acids of Human claudin 14
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	26 kDa
Gene Name:	claudin 14
Database Link:	NP_036262 Entrez Gene 56173 Mouse Entrez Gene 23562 Human O95500



[View online »](#)

Background:

Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets; forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet; with complementary grooves in the inwardly facing extracytoplasmic leaflet. The protein encoded by this gene; a member of the claudin family; is an integral membrane protein and a component of tight junction strands. The encoded protein also binds specifically to the WW domain of Yes-associated protein. Defects in this gene are the cause of an autosomal recessive form of nonsyndromic sensorineural deafness. It is also reported that four synonymous variants in this gene are associated with kidney stones and reduced bone mineral density. Several transcript variants encoding the same protein have been found for this gene.

Synonyms:


DFNB29

Protein Families:

Druggable Genome, Transmembrane

Protein Pathways:

Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction

Product images:

Gel: 10%SDS-PAGE

Lysate: 30 µg

Lane: Mouse pancreas tissue

Primary antibody: [TA322417] (CLDN14 Antibody)
at dilution 1/200Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution

Exposure time: 50 seconds