

Product datasheet for **TA322630S**

Alanine Transaminase (GPT) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 10-50 Positive control: Human prostate cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 1-14 amino acids of Human glutamic-pyruvate transaminase (alanine aminotransferase)
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54.6 kDa
Gene Name:	glutamic-pyruvate transaminase (alanine aminotransferase)
Database Link:	NP_005300 Entrez Gene 2875 Human P24298
Background:	This gene encodes cytosolic alanine aminotransaminase 1 (ALT1); also known as glutamate-pyruvate transaminase 1. This enzyme catalyzes the reversible transamination between alanine and 2-oxoglutarate to generate pyruvate and glutamate and; therefore; plays a key role in the intermediary metabolism of glucose and amino acids. Serum activity levels of this enzyme are routinely used as a biomarker of liver injury caused by drug toxicity; infection; alcohol; and steatosis. A related gene on chromosome 16 encodes a putative mitochondrial alanine aminotransaminase.




[View online »](#)

Synonyms: AAT1; ALT1; GPT1


Protein Families: Druggable Genome

Protein Pathways: Alanine, aspartate and glutamate metabolism, Metabolic pathways

Product images:



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA322630] (GPT Antibody) at dilution 1/30 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using [TA322630] (GPT Antibody) at dilution 1/30, treated with synthetic peptide. (Original magnification: $\times 200$)