

Product datasheet for **TA322705**

PERK (EIF2AK3) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 100-200 Positive control: Human colorectal cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to C terminal 250 amino acids of human eukaryotic translation initiation factor 2-alpha kinase 3
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	eukaryotic translation initiation factor 2 alpha kinase 3
Database Link:	NP_004827 Entrez Gene 13666 Mouse Entrez Gene 29702 Rat Entrez Gene 9451 Human Q9NZJ5
Background:	The protein encoded by this gene phosphorylates the alpha subunit of eukaryotic translation-initiation factor 2 (EIF2); leading to its inactivation; and thus to a rapid reduction of translational initiation and repression of global protein synthesis. It is a type I membrane protein located in the endoplasmic reticulum (ER); where it is induced by ER stress caused by malformed proteins. Mutations in this gene are associated with Wolcott-Rallison syndrome.
Synonyms:	PEK; PERK; WRS
Protein Families:	Druggable Genome, Protein Kinase, Secreted Protein, Transmembrane




[View online »](#)

Protein Pathways: Alzheimer's disease


Product images:




Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using TA322705 (EIF2AK3 Antibody) at dilution 1/90 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using TA322705 (EIF2AK3 Antibody) at dilution 1/90, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA322705 (EIF2AK3 Antibody) at dilution 1/90 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA322705 (EIF2AK3 Antibody) at dilution 1/90, treated with fusion protein. (Original magnification: $\times 200$)