

Product datasheet for **TA322867S**

TXN Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: A172 and hela cells IHC: 25-100 Positive control: Human cervical cancer Predicted cell location: Cytoplasm, Nucleus
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 91-105 amino acids of human thioredoxin
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	12 kDa
Gene Name:	thioredoxin
Database Link:	NP_003320 Entrez Gene 7295 Human P10599



[View online »](#)

Background:


The protein encoded by this gene acts as a homodimer and is involved in many redox reactions. The encoded protein is active in the reversible S-nitrosylation of cysteines in certain proteins; which is part of the response to intracellular nitric oxide. This protein is found in the cytoplasm. Two transcript variants encoding different isoforms have been found for this gene. Participates in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyzes dithiol-disulfide exchange reactions.

Synonyms:

TRDX; TRX; TRX1

Protein Families:

Druggable Genome

Product images:

Gel: 12%SDS-PAGE


Lysate: 40 μ g

Lane 1-2: A172 cells
hela cells


Primary antibody: [TA322867] (TXN Antibody) at
dilution 1/300

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution


Exposure time: 10 seconds




Immunohistochemistry of paraffin-embedded
Human cervical cancer tissue using [TA322867]
(TXN Antibody) at dilution 1/30 (Original
magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA322867] (TXN Antibody) at dilution 1/30, treated with synthetic peptide. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using [TA322867] (TXN Antibody) at dilution 1/30 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using [TA322867] (TXN Antibody) at dilution 1/30, treated with synthetic peptide. (Original magnification: $\times 200$)