

Product datasheet for TA323150

IGFBP7 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 50-200

Positive control: Human colon cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Synthetic peptide corresponding to a region derived from 270-282 amino acids of Human

insulin-like growth factor binding protein 7

Formulation: PBS pH7.3, 0.05% NaN3, 50% glycerol

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: insulin like growth factor binding protein 7

Database Link: NP 001544

Entrez Gene 29817 MouseEntrez Gene 3490 Human

Q16270

Background: This gene encodes a member of the insulin-like growth factor (IGF)-binding protein (IGFBP)

family. IGFBPs bind IGFs with high affinity; and regulate IGF availability in body fluids and tissues and modulate IGF binding to its receptors. This protein binds IGF-I and IGF-II with relatively low affinity; and belongs to a subfamily of low-affinity IGFBPs. It also stimulates prostacyclin production and cell adhesion. Alternatively spliced transcript variants encoding different isoforms have been described for this gene; and one variant has been associated

with retinal arterial macroaneurysm.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Synonyms: AGM; FSTL2; IBP-7; IGFBP-7v; IGFBPRP1; MAC25; PSF; RAMSVPS; TAF

Protein Families: Secreted Protein

Product images:




Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA323150 (IGFBP7 Antibody) at dilution 1/50 (Original magnification: ×200)




Immunohistochemistry of paraffin-embedded Human colon cancer tissue using TA323150 (IGFBP7 Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA323150 (IGFBP7 Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA323150 (IGFBP7 Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: ×200)