

Product datasheet for TA323354

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Sex Hormone Binding Globulin (SHBG) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 25-100

Positive control: Human cervical cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Synthetic peptide corresponding to a region derived from 68-81 amino acids of Human sex

hormone-binding globulin

Formulation: PBS pH7.3, 0.05% NaN3, 50% glycerol

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: sex hormone binding globulin

Database Link: NP 001031

Entrez Gene 20415 MouseEntrez Gene 24775 RatEntrez Gene 6462 Human

P04278

Background: This gene encodes a steroid binding protein that was first described as a plasma protein

secreted by the liver but is now thought to participate in the regulation of steroid responses. The encoded protein binds each steroid molecule as a dimer formed from identical or nearly identical monomers. The use of alternate promoters and alternatively spliced transcripts have been described. Multiple transcript variants encoding different isoforms have been

found for this gene.


Synonyms: ABP; SBP; TEBG






Protein Families: Druggable Genome, Secreted Protein

Product images:




Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA323354 (SHBG Antibody) at dilution 1/20 (Original magnification: ×200)




Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA323354 (SHBG Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA323354 (SHBG Antibody) at dilution 1/20 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using TA323354 (SHBG Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: ×200)