

## Product datasheet for **TA323461S**


### MAGEA4 Rabbit Polyclonal Antibody

#### Product data:


|                       |  |
|-----------------------|--|
| Product Type:         | Primary Antibodies   |
| Applications:         | IHC  |
| Recommended Dilution: | IHC: 10-50<br>Positive control: Human tonsil<br>Predicted cell location: Cytoplasm   |
| Reactivity:           | Human  |
| Host:                 | Rabbit   |
| Isotype:              | IgG  |
| Clonality:            | Polyclonal   |
| Immunogen:            | Full length fusion protein   |
| Formulation:          | PBS pH7.3, 0.05% NaN <sub>3</sub> , 50% glycerol   |
| Purification:         | Antigen affinity purification  |
| Conjugation:          | Unconjugated   |
| Storage:              | Store at -20°C as received.  |
| Stability:            | Stable for 12 months from date of receipt.   |
| Gene Name:            | MAGE family member A4  |
| Database Link:        | <a href="#">NP_002353</a><br><a href="#">Entrez Gene 4103 Human</a><br><a href="#">P43358</a>  |
| Background:           | This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability; suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders; such as dyskeratosis congenita. At least four variants encoding the same protein have been found for this gene. |
| Synonyms:             | CT1.4; MAGE-41; MAGE-X2; MAGE4; MAGE4A; MAGE4B   |



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA323461] (MAGEA4 Antibody) at dilution 1/10 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA323461] (MAGEA4 Antibody) at dilution 1/10, treated with fusion protein. (Original magnification: ×200)