

## Product datasheet for **TA323520**

### Prolactin (PRL) Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to a region derived from 29-227 amino acids of human prolactin
Formulation:	PBS pH7.3, 0.05% NaN <sub>3</sub> , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	prolactin
Database Link:	<a href="#">NP_000939</a> <a href="#">Entrez Gene 5617 Human P01236</a>
Background:	This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues; including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis; and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein.
Synonyms:	GHA1
Protein Families:	Druggable Genome, Secreted Protein




[View online »](#)

**Protein Pathways:** Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-receptor interaction

**Product images:**



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA323520 (PRL Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA323520 (PRL Antibody) at dilution 1/30, treated with fusion protein. (Original magnification: ×200)