

Product datasheet for **TA323584**

ALDH9A1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: Human fetal liver tissue and hela cells, Human fetal brain tissue IHC: 50-200 Positive control: Human breast cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 18-32 amino acids of human aldehyde dehydrogenase 9 family, member A1
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54 kDa
Gene Name:	aldehyde dehydrogenase 9 family member A1
Database Link:	NP_000687 Entrez Gene 56752 Mouse Entrez Gene 64040 Rat Entrez Gene 223 Human P49189



[View online »](#)

Background:

This protein belongs to the aldehyde dehydrogenase family of proteins. It has a high activity for oxidation of gamma-aminobutyraldehyde and other amino aldehydes. The enzyme catalyzes the dehydrogenation of gamma-aminobutyraldehyde to gamma-aminobutyric acid (GABA). This isozyme is a tetramer of identical 54-kD subunits. Converts gamma-trimethylaminobutyraldehyde into gamma-butyrobetaine. Catalyzes the irreversible oxidation of a broad range of aldehydes to the corresponding acids in an NAD-dependent reaction.

Synonyms:


ALDH4; ALDH7; ALDH9; E3; TMABADH

Protein Families:

Druggable Genome

Protein Pathways:

Arginine and proline metabolism, Ascorbate and aldarate metabolism, beta-Alanine metabolism, Butanoate metabolism, Fatty acid metabolism, Glycerolipid metabolism, Glycolysis / Gluconeogenesis, Histidine metabolism, Limonene and pinene degradation, Lysine degradation, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism, Tryptophan metabolism, Valine, leucine and isoleucine degradation

Product images:

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane 1-3: Human fetal liver tissue
hela cells

Human fetal brain tissue


Primary antibody: TA323584 (ALDH9A1 Antibody)
at dilution 1/240

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution


Exposure time: 30 seconds




Immunohistochemistry of paraffin-embedded Human breast cancer tissue using TA323584 (ALDH9A1 Antibody) at dilution 1/50 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using TA323584 (ALDH9A1 Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA323584 (ALDH9A1 Antibody) at dilution 1/50 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA323584 (ALDH9A1 Antibody) at dilution 1/50, treated with synthetic peptide. (Original magnification: $\times 200$)