

## Product datasheet for **TA323587**

### alpha 1b Adrenergic Receptor (ADRA1B) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 100-500 WB positive control: Human lung cancer tissue
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 500-515 amino acids of Human Alpha-1B adrenergic receptor
Formulation:	PBS pH7.3, 0.05% NaN <sub>3</sub> , 50% glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57 kDa
Gene Name:	adrenoceptor alpha 1B
Database Link:	<a href="#">NP_000670</a> <a href="#">Entrez Gene 11548 Mouse</a> <a href="#">Entrez Gene 24173 Rat</a> <a href="#">Entrez Gene 147 Human</a> <a href="#">P35368</a>



[View online »](#)

**Background:**

Alpha-1-adrenergic receptors (alpha-1-ARs) are members of the G protein-coupled receptor superfamily. They activate mitogenic responses and regulate growth and proliferation of many cells. There are 3 alpha-1-AR subtypes: alpha-1A; -1B and -1D; all of which signal through the Gq/11 family of G-proteins and different subtypes show different patterns of activation. This gene encodes alpha-1B-adrenergic receptor; which induces neoplastic transformation when transfected into NIH 3T3 fibroblasts and other cell lines. Thus; this normal cellular gene is identified as a protooncogene. This gene comprises 2 exons and a single large intron of at least 20 kb that interrupts the coding region.

**Synonyms:**

ADRA1; ALPHA1BAR

**Protein Families:**

Druggable Genome, GPCR, Transmembrane

**Protein Pathways:**

Calcium signaling pathway, Neuroactive ligand-receptor interaction, Vascular smooth muscle contraction

**Product images:**

Gel: 10%SDS-PAGE

Lysate: 40 µg

Lane: Human lung cancer tissue

Primary antibody: TA323587 (ADRA1B Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute