

## Product datasheet for **TA323902S**


### ADCK4 (COQ8B) Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human colon cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to a region derived from 192-423 amino acids of human aarF domain containing kinase 4
Formulation:	PBS pH7.3, 0.05% NaN <sub>3</sub> , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	aarF domain containing kinase 4
Database Link:	<a href="#">NP_079152</a> <a href="#">Entrez Gene 76889 Mouse</a> <a href="#">Entrez Gene 308453 Rat</a> <a href="#">Entrez Gene 79934 Human</a> <a href="#">Q96D53</a>
Background:	This gene encodes a protein with two copies of a domain found in protein kinases. The encoded protein has a complete protein kinase catalytic domain, and a truncated domain that contains only the active and binding sites of the protein kinase domain, however, it is not known whether the protein has any kinase activity. Multiple transcript variants encoding different isoforms have been found for this gene.?
Synonyms:	NPHS9
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA323902] (COQ8B Antibody) at dilution 1/25 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using [TA323902] (COQ8B Antibody) at dilution 1/25, treated with fusion protein. (Original magnification:  $\times 200$ )