

Product datasheet for **TA324102S**

MMP11 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human colon and ovarian cancer tissue
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a region derived from 470-488 amino acids of Human matrix metalloproteinase 11 (stromelysin 3)
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55 kDa
Gene Name:	matrix metalloproteinase 11
Database Link:	NP_005931 Entrez Gene 17385 Mouse Entrez Gene 4320 Human P24347
Background:	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes; such as embryonic development; reproduction; and tissue remodeling; as well as in disease processes; such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However; the enzyme encoded by this gene is activated intracellularly by furin within the constitutive secretory pathway. Also in contrast to other MMP's; this enzyme cleaves alpha 1-proteinase inhibitor but weakly degrades structural proteins of the extracellular matrix.




[View online »](#)

Synonyms: SL-3; ST3; STMY3

Protein Families: Druggable Genome, Protease

Product images:



Gel: 10%SDS-PAGE

Lysate: 60 μ g

Lane 1-2: Human colon tissue

Human ovarian cancer tissue

Primary antibody: [TA324102] (MMP11 Antibody)
at dilution 1/600

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution

Exposure time: 20 seconds