

Product datasheet for **TA324263S**

MCM3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: K562, Jurkat, hela and 293T cells IHC: 25-100 Positive control: Human lung cancer Predicted cell location: Cytoplasm, Nucleus
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein corresponding to C terminal 300 amino acids of human minichromosome maintenance complex component 3
Formulation:	PBS pH7.3, 0.05% NaN ₃ , 50% glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	91 kDa
Gene Name:	minichromosome maintenance complex component 3
Database Link:	NP_001257401 Entrez Gene 17215 Mouse Entrez Gene 4172 Human P25205



[View online »](#)

Background:

The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are involved in the initiation of eukaryotic genome replication. The hexameric protein complex formed by MCM proteins is a key component of the pre-replication complex (pre_RC) and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. This protein is a subunit of the protein complex that consists of MCM2-7. It has been shown to interact directly with MCM5/CDC46. This protein also interacts with and is acetylated by MCM3AP; a chromatin-associated acetyltransferase. The acetylation of this protein inhibits the initiation of DNA replication and cell cycle progression. Two transcript variants encoding different isoforms have been found for this gene.

Synonyms:


HCC5; P1-MCM3; P1.h; RLFB

Protein Families:

Druggable Genome, Stem cell - Pluripotency, Transcription Factors

Protein Pathways:

Cell cycle, DNA replication

Product images:

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane 1-4: K562 cells

Jurkat cells


hela cells

293T cells


Primary antibody: [TA324263] (MCM3 Antibody)
at dilution 1/550

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution

Exposure time: 1 minute



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using [TA324263] (MCM3 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using [TA324263] (MCM3 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: ×200)