

Product datasheet for TA324323

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

TUBB1 Mouse Monoclonal Antibody [Clone ID: 86CT4.4.4]

Product data:

Product Type: Primary Antibodies

Clone Name: 86CT4.4.4

Applications: IF, IHC, WB

Recommended Dilution: WB: 1:100~500, IHC: 1:10~50, IF: 1:10~50

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: This monoclonal antibody is generated from mice immunized with a recombinant protein

human TBB1.

Formulation: Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.

Concentration: lot specific

Purification: This antibody is purified through a protein G column, followed by dialysis against PBS.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 50327 Da

Gene Name: tubulin beta 1 class VI

Database Link: NP 110400

Entrez Gene 545486 MouseEntrez Gene 81027 Human

Q9H4B7


Synonyms: beta 1; beta tubulin 1; class VI; dJ543J19.4; OTTHUMP00000031411; tubulin

Protein Families: Druggable Genome


Protein Pathways: Gap junction, Pathogenic Escherichia coli infection




Product images:



Western blot analysis of anti-TBB1 Antibody (Cat. #TA324323) in mouse brain tissue, Y79, CEM and 293 lysates (35ug/lane). TBB1 (arrow) was detected using the purified Mab. (1ug/ml)



Formalin-fixed and paraffin-embedded human brain reacted with TBB1 Antibody, which was peroxidase-conjugated to the secondary antibody, followed by DAB staining. This data demonstrates the use of this antibody for immunohistochemistry; clinical relevance has not been evaluated.



Confocal immunofluorescent analysis of TBB1 Antibody (Cat#TA324323) with HepG2 cell followed by Alexa Fluor® 488-conjugated goat anti-mouse lgG (green). DAPI was used to stain the cell nuclear (blue).