


## Product datasheet for **TA324508S**

### LZTS1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This LZTS1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 77-105 amino acids from the N-terminal region of human LZTS1.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	66613 Da
Gene Name:	leucine zipper, putative tumor suppressor 1
Database Link:	<a href="#">NP_066300</a> <a href="#">Entrez Gene 211134 Mouse</a> <a href="#">Entrez Gene 11178 Human</a> <a href="#">Q9Y250</a>
Synonyms:	F37; FEZ1
Protein Families:	Druggable Genome

[View online »](#)

**Product images:**


LZTS1 Antibody (N-term) (Cat. #[TA324508])  
western blot analysis in A2058 cell line and  
mouse brain tissue lysates (35ug/lane). This  
demonstrates the LZTS1 antibody detected the  
LZTS1 protein (arrow).