


Product datasheet for **TA324596S**

ZNF90 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This ZNF90 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 84-112 amino acids from the Central region of human ZNF90.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	69058 Da
Gene Name:	zinc finger protein 90
Database Link:	NP_009069 Entrez Gene 7643 Human Q03938
Synonyms:	HTF9

[View online »](#)

Product images:


ZNF90 Antibody (Center) (Cat. #[TA324596])
western blot analysis in CEM cell line lysates
(35ug/lane). This demonstrates the ZNF90
antibody detected the ZNF90 protein (arrow).